

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\DTI_Gradient_Test\dti_nopf_x2_ortho
 TA:2:16 PAT:2 Voxel size:3.0x3.0x3.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5600 ms
TE	79.0 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	72
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Fat sat. mode	Weak
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	108 mm
Frequency 1H	123.255749 MHz
Correction factor	1
AddCSaCSatNS 1H	66.311 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1336 Hz/Px
Optimization	None
Free echo spacing	Off
Echo spacing	0.93 ms
EPI factor	72
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\DTI_Gradient_Test\dti_nopf_x2_roll
 TA:2:16 PAT:2 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5600 ms
TE	79.0 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	72
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Fat sat. mode	Weak
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	108 mm
Frequency 1H	123.255749 MHz
Correction factor	1
AddCSaCSatNS 1H	66.311 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1336 Hz/Px
Optimization	None
Free echo spacing	Off
Echo spacing	0.93 ms
EPI factor	72
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\DTI_Gradient_Test\dti_nopf_x2_pitch
 TA:2:16 PAT:2 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5600 ms
TE	79.0 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	72
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Fat sat. mode	Weak
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	108 mm
Frequency 1H	123.255749 MHz
Correction factor	1
AddCSaCSatNS 1H	66.311 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1336 Hz/Px
Optimization	None
Free echo spacing	Off
Echo spacing	0.93 ms
EPI factor	72
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\DTI_Gradient_Test\dti_nopf_x2_yaw
 TA:2:16 PAT:2 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5600 ms
TE	79.0 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	72
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Fat sat. mode	Weak
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	108 mm
Frequency 1H	123.255749 MHz
Correction factor	1
AddCSaCSatNS 1H	66.311 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1336 Hz/Px
Optimization	None
Free echo spacing	Off
Echo spacing	0.93 ms
EPI factor	72
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\DTI_Gradient_Test\dti_nopf_x2_2axis
 TA:2:16 PAT:2 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	5600 ms
TE	79.0 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	72
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Fat sat. mode	Weak
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	108 mm
Frequency 1H	123.255749 MHz
Correction factor	1
AddCSaCSatNS 1H	66.311 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1336 Hz/Px
Optimization	None
Free echo spacing	Off
Echo spacing	0.93 ms
EPI factor	72
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off